

# TENSIMETER AIR RAKSA FUJITO



## Spesifikasi

Versions : Aneroid, Desk, Wall, Stand and anaesthetic models

Type of indication : Round column

Scale graduation : Increments of 2 mm Hg

Range shown on scale : 0 to 300 mm Hg

Measurement range : 0 to 300 mm Hg

Pointer movement: No stop pin

Pressure Generation: Bulb

Pressure reduction: Air- release valve that can be regulated

**PT. CENTRAL BANDUNG RAYA**

Jl. Kom. Ud. Supadio No. 31 Bandung 40174

Telp. 022-6014183, 6071214, 6077413, 6032021 Fax. 022-6032021

Email : [cbr@bdg.centrin.net.id](mailto:cbr@bdg.centrin.net.id)

Blog : [www.centralbandungraya.wordpress.com](http://www.centralbandungraya.wordpress.com)

Website : [www.centralbandungraya.com](http://www.centralbandungraya.com)